Table 1

| Light intensity (LED-/halogen lamp) | Polymerisation time (per layer) |
|-------------------------------------|---|
| ≥ 300 mW/cm² (LED-lamp) | 20 s |
| ≥ 500 mW/cm² (halogen lamp) | 40 s A1, A2, A3, A3.5, BL A4, OA2 |
| ≥ 1000 mW/cm² | 10 s A1, A2, A3, BL A3.5, A4, OA2 |

INSTRUCTIONS FOR USE

EN

MD EU Medical device

In accordance with DIN EN ISO 4049

Product description

Visalys® Flow is a light-curing flowable restorative combining proven composite technology with innovative nano-technology. The radiopaque Visalys® Flow contains 80% w/w inorganic fillers ($\Delta 65.6\%$ vol.) in a methacrylate matrix and cures under blue light (halogen/LED). It is suited for restorative therapy in anterior and posterior areas. Because of its low viscosity, the material exhibits very good wetting and adaptation to cavity walls. Visalys® Flow is used with a dentine enamel bond.

Shades

A1, A2, A3, A3.5, A4, OA2, BL

Indications

- Filling minimally invasive cavities of all classes
- Filling small class I cavities and extended fissure sealing
- Filling class II-V cavities including V-shaped defects and cervical caries
- Blocking out undercuts
- Lining or coating cavities
- Repairing fillings and veneers
- Luting translucent prosthetic pieces (e.g. full ceramic crowns, etc.)

Contraindications

Visalys® Flow contains methacrylate and BHT. Visalys® Flow should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group

Visalys® Flow is suitable for use in all patients without any age or gender restrictions.

Performance features

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User

Visalys® Flow should only be used by a professionally trained dental practitioner.

Use

Visalys® Flow is a material which adapts very well to the cavity, helps to avoid air enclosures and accelerates restorations.

Preparation

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection

Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

Cavity preparation

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. Bevel enamel margins. Observe special preparation forms in extended fissure sealing and deciduous teeth! Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Creating a dry working field

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement

Take appropriate measures for areas where the applied Visalys® Flow can flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

Pulp protection

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding materials

Visalys® Flow is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application

Lock the application tip onto the syringe with a clockwise turn and apply Visalys® Flow directly. For fillings of more than 2 mm, apply and polymerise in layers. In order to ensure the function of the non-dripping syringe do not pull back the plunger during or after use.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing

Conventional polymerisation devices are suitable for light-curing the material. The polymerisation time per layer (Polymerisation time) as a function of the light intensity (Light intensity) can be found in Table 1 (Table 1).

Hold the tip of the curing light as close as possible to the surface of the filling. If the distance is more than 5mm, the curing depth will be reduced. Insufficient curing can lead to discolouration and pulp irritation.

Finishing

The filling can be finished and polished immediately after removing the matrices (e.g. fine and extra-fine grain diamonds, polishing discs).

During finishing and polishing, the filling must be cooled. The margin of the filling, or even better, the entire tooth should be fluoridated as a final step.

Warnings, precautionary measures

- No known side effects. Hypersensitive patients may experience sensitivity.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of Visalys® Flow. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with Visalys® Flow.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content)

Strontrium aluminum borosilicate glass, Barium aluminium borosilicate glass, silicon dioxide, HEDMA, BisGMA, TEGDMA, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method

Storage at 4°C–28°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal

Dispose of the product in accordance with local regulations.

Reporting obligation

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of Visalys® Flow must be reported to the manufacturer and the responsible authority.

Note

The Summary of Safety and Clinical Performance of Visalys® Flow can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

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