1.3 Superfici adesive in composto

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Englisch

Product description / Indication

Visalys® Restorative Primer with adhesive monomers (10-MDP, silane methacrylate) is for building up adhesive surfaces on restoration materials (metal, composite, as well as silicate and oxide ceramics). Visalys® Restorative Primer is a single-component primer. Do not light cure Visalys® Restorative Primer. BPA or BPA precursors are not used in the manufacturing process of this device.

Indications for use

Restoration primer for treatment of metal/silicate ceramics, oxide ceramic and composite surfaces

- to cement prosthetic restorations
- to repair fractured restorations with composite materials

Contraindications

In case of known allergies to methacrylate monomers or other constituents of Visalys® Restorative Primer.

If there are no adequate means of ensuring drainage.

Warnings

Allergic reactions to the product cannot be ruled out in sensitive persons. Sensitization through contact with skin / mucous tissue is possible.

The use of hydrogen peroxide solutions, materials containing eugenol or retraction paste containing iron(III) sulfate can impair adhesion.

Keep out of reach for children.

Cautions

Do not use after the expiration date.

Do not swallow material. If the material is inadvertently swallowed: Consult a medical doctor if problems arise or persist. Take the appropriate protective measures (e.g. wear gloves or safety glasses). Avoid contact of skin and mucous tissue with the material.

Avoid contact with the eyes. If there is accidental contact with the eyes: rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Avoid contact with clothing, as permanent damage to clothes cannot be excluded.

Safety data sheets are available at www.kettenbach.com.

General notes

Given refrigerator storage, the product has to have reached room temperature prior to use.

Visalys® Restorative Primer is intended for surface treatment of prosthetic restorations, indirect core build-ups and implant abutments and not for pretreatment of natural tooth hard substance.

Visalys® Restorative Primer should be processed quickly, as the constituents can evaporate.

The adhesive surface must be clean and not contaminated, e.g. with saliva, blood or water, prior to applying Visalys® Restorative Primer.

The surfaces of some types of posts can be attacked by sandblasting or hydrofluoric acid. The specifications of the root post manufacturer on surface treatment must be observed.

For cementation on an implant abutment or on an indirect core build-up, the specifications of the respective material manufacturer must always be considered. Unless prescribed otherwise in the relevant material manufacturer’s instructions for use, pre-treat the entire adhesive surfaces of the abutment / core build-up dependent on the material, as described under "Pre-treatment of the restoration" in the Visalys® CemCore instructions for use.

The emergence profiles and screw channels, as applicable, must be protected by suitable measures (e.g. wax, silicone covers).

The cementation composite (e.g. Visalys® CemCore) must then be applied according to the instructions for use of the product.

Use as primer on restoration materials in connection with cementation composites

After the try-in, the surfaces of the restoration contaminated e.g. with saliva, blood or silicone residues and any Visalys® CemCore Try In Paste, have to be cleaned with distilled water spray or in the ultrasonic bath and then dried.

All the surfaces must be free of contamination. In the case of restorations made of non-precious metal or oxide ceramics, there

must be no cleaning with phosphoric acid.

Unless prescribed otherwise in the instructions for use for the restoration material, root post or material for the indirect core build-up, the following procedure should be selected for pretreatment:

- Adhesive surfaces of noble metal, metal and oxide ceramics (zirconium oxide and aluminum oxide ceramics)
 - Sand blasting of the surfaces with aluminum oxide powder (≤ 50 µm). Observe the restoration material manufacturer’s blasting parameters.
 - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.
 - DO NOT CLEAN WITH PHOSPHORIC ACID.**
 - Dry with oil-free air.
- Adhesive surfaces of silicate ceramic (feldspar ceramic, glass ceramics) and hybrid ceramics
 - Etch the surfaces with hydrofluoric acid according to the specifications in the instructions for use of the restoration material.
 - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.
 - Dry with oil-free air.
- Adhesive surfaces of composite
 - Sand blasting of the surfaces with aluminum oxide powder (≤ 50 µm). Observe the restoration material manufacturer’s blasting parameters.
 - Then clean thoroughly with alcohol.
 - Dry with oil-free air.

- Drop Visalys® Restorative Primer into a preparation tray. Light protection is not necessary. Close primer bottles carefully after each use.
- Apply Visalys® Restorative Primer with a disposable applicator or brush onto the adhesive surfaces and allow it to work in for 60 sec. Then blow out and dry with a gentle jet of air. Any excess can be removed with a cotton pellet. The adhesive bond to precious metal alloys can be enhanced by using a special precious metal primer. Visalys® Restorative Primer should be processed quickly, as the constituents can evaporate.

Use as primer to repair fractured restorations

- Roughen the entire fractured surfaces of the restoration material by blasting or using a suitable rotary instrument. Observe the restoration material manufacturer’s blasting parameters and other measures, as applicable.
- Clean and dry the fractured surface. Ensure a clean and dry operating field (rubber dam) during the entire application.
- Drop Visalys® Restorative Primer into a preparation tray. Light protection is not necessary.
- Apply Visalys® Restorative Primer with a disposable applicator or brush onto the adhesive surfaces and allow it to work in for 60 sec. and then blow out and dry with a gentle oil-free jet of air. Any excess can be removed with a cotton pellet. The adhesive bond to precious metal alloys can be enhanced by using a special precious metal primer. Visalys® Restorative Primer should be processed quickly, as the constituents can evaporate.
- Carry out the repair with a composite (e.g. Visalys® CemCore) according to the manufacturer’s specifications.

Note

If part of the fractured surface consists of natural tooth substance, they have to be pretreated with adhesive. The specifications of the adhesive manufacturer must be observed.

Cleaning and disinfection

The Visalys® Restorative Primer bottle can be cleaned with wipe disinfection using non-oxidizing disinfectants.

Disinfection of the bottle in accordance with national guidelines.

Observe standard practice hygiene measures.

Storage Conditions

Store dry, protected from sunlight. The product has to be kept in the refrigerator at 2-8°C/35-46°F before first use.

After the first use, Visalys® Restorative Primer has to be used within 6 months. Avoid extreme temperatures. Given refrigerator storage, the product has to have reached room temperature prior to use.

Do not use after the expiration date.

Constituents

The main components of the Visalys® Restorative Primer are ethanol, acidic adhesive monomer (10-MDP) and silane methacrylate.

Disposal

Dispose of the product in accordance with local regulations.

Trademarks

Visalys® Restorative Primer, Visalys® CemCore Try In Paste and Visalys® CemCore are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Not all Kettenbach products are available in all countries. For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

Limitation of liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal basis asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to

technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use of the product.

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